- Make sure that the use of the device is not unpleasant and does not impair dental health. If in doubt, contact a dentist. If you have dental problems, wear braces, have dental implants or crowns, consult your dentist or physician before using the product.
- Do not open the product. Any warranty is void if the product is opened. Improper repair is dangerous for the user.
- Never use a product that has been opened or damaged.
- Do not wash the product in the dishwasher.
- People with pacemakers or other electronic implants should consult their physician beforehand.
- The product should not be used for commercial purposes.
- Keep the packaging out of the reach of children.
- The manufacturer is not liable for any failure to observe these operating instructions, for repairs that are not carried out in accordance with the state of the art, or for improper use of the device.

Attention.

Please note that the warranty does not apply if water is found in the handle or if the handle has not been used according to these instructions. It is therefore important to close the connector cover and not to leave the handle within reach of a sink, shower or bath.

Note for USA:

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures :

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Note for Canada : CAN ICES-3 (B) / NMB-3(B)

This device complies with Canadian ICES-003 Class B and RSS-210.

Operation is subject to the following two conditions :

- (1) This device may not cause interference, and
- (2) This device must accept any interference, including interference that may cause undesired operation of the device.

Appliance

European directive 2012/19/EU on Waste from Electrical and Electronic Equipment (WEEE), requires that used household appliances are not thrown into the normal municipal waste stream. Used appliances must be collected separately in order to optimize the rate of recovery and recycling of materials that compose them, and to reduce the impact on human health and on the environment. The crossed bin symbol is affixed to all the products to remind you of the obligations of separated collection

Battery

Do not dispose of used batteries with household waste. They should be disposed of at a collection point. You can obtain more information from your local waste management company or your local council.

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Site : www.y-brush.com

Operating instructions

Please read the operating instructions carefully before using the device for the first time. Toothbrushes can be used by children and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved. Cleaning and user maintenance shall not be made by children without supervision. The device should be used only with the power supply or USB charger provided with the product and the brush sold by Y-Brush. Connect the charger only to an AC outlet (100-240V). The product should only be operated with the safety extra-low voltage corresponding to the marking on the product. Children shall not play with the appliance. This appliance contains battery that is non-replaceable.

IMPORTANT SAFEGUARDS

When using electrical products, especially when children are present, basic safety precautions should always be followed, including the following:

READ ALL INSTRUCTIONS BEFORE USING

DANGER- To reduce the risk of electrocution:

1. Do not use while bathing.
2. Do not place or store product where it can fall or be pulled into a tub or sink.
3. Do not place in or drop into water or other liquid.
4. Do not reach for a product that has fallen into water. Unplug immediately.

WARNING

To reduce the risk of burns, electrocution, fire, or injury to persons :

1. Close supervision is necessary when this product is used by, on, or near children or invalids.
2. Use this product only for its intended use as described in this manual. Do not use attachments not recommended by the manufacturer.
3. Never operate this product if it has a damaged cord or plug, if it is not working properly, if it has been dropped or damaged, or dropped into water. Return the product to a service center for examination and repair.
4. Keep the cord away from heated surfaces.
5. Never use while sleeping or drowsy.
6. Do not use outdoors or operate where aerosol (spray) products are being used or where oxygen is being administered.

INTRODUCTION

Y-Brush is a device designed and developed by oral health specialists to facilitate dental brushing for both children and adults. It allows you to brush your teeth in a simple, effective way, and in a few seconds thanks to a patented system. This device is composed of a flexible brush that fits perfectly in the jaw. Thanks to the nylon fibers lining the brush and the vibrations of the handle, the teeth are brushed simultaneously and efficiently. Remember that the recommendations in terms of oral health are to brush your teeth at least twice a day, in the morning and in the evening. An effective brushing helps to prevent the consequences related to a bad oral health, (periodontal, cardiovascular diseases, etc.) If helps to remove plaque, but is not intended to remove tartar deposits. It is recommended to see a dental surgeon every six months.

DIRECTION FOR USE

+ CHARGING DEVICE

- Before each use of the Y-Brush, make sure it is fully charged. In normal use, the indicator lights on the handle are lit in white. They will light in red when the battery level is too low and in green when the Y-Brush is plugged into the power supply.
- If the battery level of the handle is too low:
 - ◊ Connect the charger to a 100-240V power outlet or use the USB charger
 - ◊ Plug the charging cable into the port provided for this purpose on the handle. An illuminated halo around the on/off button lights up cyclically and indicates that the device is being charged. If the halo stays on green continuously, it indicates that the handle is fully charged. A period of 8 hours is sufficient to fully charge the handle. It will then have an autonomy of up to 3 months at a rate of 2 brushes per day.
 - ◊ The handle is ready to be used. It is imperative to close the cover that protects the connection of the recharging cable on the handle side to guarantee its watertightness.
 - ◊ The handle cannot be used during charging.

WARNING

THE BRUSH

- The appliance and its brush have the sole function of brushing teeth, they should not be used for anything else. Do not attempt to swallow all or part of the brush.
- The handle must not be inserted in the mouth without a brush.
- Make sure that the nylon fibers of the brush reach all areas of all teeth once in the mouth. If any area are not reached, brush them separately.
- Be aware of the risk of pinching when the brush is attached to the handle.
- The brushes are individual and have a maximum service life of 6 months under normal use. If there is any doubt as to whether the brush belongs to a third party and/or can be used by a third party, change the brush to avoid the risk of cross-contamination.
- Replace the brush after severe impact (even if there is no visible sign of damage).

+ POWER SUPPLY & CONNECTION

- Use the unit at room temperature (0-40°C). Avoid prolonged exposure of the product to sunlight or any other heat source. Do not use the charger near water (bathtub, washbasin, etc.) or outdoors. Never immerse or rinse the charger.
- Never touch the power cable with sharp or pointed objects. Never touch the plug with wet hands, and never unplug by pulling on the cable. Regularly check the cable for cracks or damage.
- Be sure to close the cover that protects the handle charging plug connection so that it remains watertight. The handle is only splash-proof.
- Do not place in a microwave oven or any other household appliance.
- After use, rinse the brush head thoroughly, shake off excess water and store it on the brush holder with the connector facing up.
- Do not dispose of the product in the household garbage. The handle has batteries and should therefore be taken to a collection point for electronic products. The charger must also be taken to such a collection point. Do not attempt to remove the batteries.

+ OTHER IMPORTANT INSTRUCTIONS

- **CAUTION:** Not suitable for children under 4 years old. The elements, in both the adult and child version, contains small parts that may cause a choking hazard. Therefore, they must not be used by a child without adult supervision.